

Certificate of Quality

PRODUCT: Sodium Chloride
Cat. No.: BR460-1KG
Lot No.:
COMPONENTS: Sodium Chloride
DESCRIPTION: BIOTECHNOLOGY GRADE
CAS NO. 7647-14-5
Formula NaCl
Molecular Weight 58.44

TEST	SPECIFICATION	DISPOSITION
Calcium	<= 0.005 %	PASS
DNase	NONE	PASS
Heavy Metals	<= 0.0005 %	PASS
Iodide	<= 0.002 %	PASS
Iron	<= 0.0002 %	PASS
Loss on Drying	<= <0.5 %	PASS
Magnesium	<= 0.005 %	PASS
Nitrogen Compounds	<= 0.001 %	PASS
pH (5%, Water) @25C	5.0 - 9.0	PASS
Phosphate	<= 0.0005 %	PASS
Potassium	<= 0.005 %	PASS
Protease	NONE	PASS
Purity	>= 99.9 %	PASS
RNase	NONE	PASS
Sulfate	<= 0.0004 %	PASS

Analysis may have been rounded to significant digits in specification limits.

Product meets analytical specifications of the grades listed.

STORAGE: Store at room temperature.

CAUTION: Avoid contact with eyes, skin and clothing.
Wash all areas of contact with copious amounts of water.
For laboratory or further manufacturing use only.