Certificate of Quality

PRODUCT: TG Buffer(10X)
Cat. No.: BR130-1L
Lot No.:
COMPONENTS: 0.25M Tris-Base
1.92M Glycine
DESCRIPTION: Results: Pass

STORAGE: Store at room temperature.

CAUTION: Avoid contact with eyes, skin and clothing.
Wash all areas of contact with copious amounts of water. For laboratory or further manufacturing use only.

2F., No. 253,
Sec. 4, Chung-Ching N. Rd., Shihlin District,
Taipei City 111, Taiwan, R.O.C.

TEL: 886-2-2811-8200
FAX: 886-2-2811-8628
rainbow@rainbowbiotech.com.tw